


STANTON

MODEL C.B.6

"University"

To meet the demand of colleges, technical institutes and works laboratories for a robust and accurate analytical balance, we now offer the "University" Balance, an inexpensive yet high grade instrument incorporating many refinements usually reserved for higher priced instruments.



***"A balance of outstanding value and utility
at a most competitive price"***

CAPACITY : 200 Grammes.

SENSITIVITY : 0.1 Milligramme.

BEAM : 12.7 cm.—graduated 100 divisions.
(Rider weight—5 mg.)

BEAM ARRESTMENT : "Fall away" type.

PANS : Chromium plated—slightly concave.

CASE : Finely polished mahogany, with side doors and counterpoised front slide.

RIDER : Single rod type ensuring, perfectly smooth action.

LEVELLING : Circular spirit level.

PAN SUPPORTS : Adjustable.

BASE : Heavy black plate glass.

RELEASE : Side action—knurled head.

KNIFE EDGES & PLANES : Selected agate.

Stanton Instruments Limited



119 OXFORD STREET - LONDON - W.1

Telephone GERRARD 7533/4

Printed in England, S. & O. London