


# STANTON

# HEAVY ★ ★ DUTY



**MODEL H.D. 2.**

*Fitted throughout with Synthetic Sapphire (Corundum) Planes*

**CAPACITY :** 5 Kilo.

**SENSITIVITY :** 20 milligrammes.

**BEAM :** 16". Specially extruded duralumin.

**BEAM ARRESTMENT :** Beam and Stirrups located by special nickel chrome points.

**PANS :** 9" diameter. Distance between wires :  $9\frac{1}{4}$ ".  
Height from pan to shoulders : 14".

**KNIFE EDGES :** Specially selected Agate.

**RELEASE :** Lever action.

*Balance mounted on polished mahogany base, complete with levelling screws and spirit level.*

**MODEL H. D. 1. :** As above but fitted with AGATE planes.

**MODEL H. D. 4. :** As model H. D. 1. Capacity : 2 Kilo. Sensitivity : 10 mg.

**MODEL H. D. 5. :** As model H. D. 4. but fitted with **SYNTHETIC SAPPHIRE (CORUNDUM)** planes.

The craftsmanship embodied in this instrument is of the very highest order, enhancing the already high reputation gained by Stanton balances throughout the world

**ALL KNIFE EDGES ARE FULLY RELIEVED WHEN BALANCE IS NOT IN ACTION**



**AGENTS**

**South Africa**

Macdonald Adams Co., P.O. Box 68, Johannesburg

**Malaya and Singapore**

Associated Inst. Mfs. Ltd., 3 Malacca St., Singapore

**Argentine and Uruguay**

London & Scandinavian Metallurgical Co., Ltd., London

**Egypt**

S.M.A. Azab & Co., 12 Antikhana Street, Cairo

## *Stanton Instruments Limited*

119 OXFORD STREET · LONDON · W.1 Telephone GERRARD 7533