


● Balances
of Precision

●
Instruments
Instanton



Metric

Imperial

Imperial

Metric



Capacity
Sensitivity

5 kg
4.0 mg

11 lb
0.06 grain

25 kg
15 mg

56 lb
0.2 grain



Scale piece
40 divisions
Knife edges
Selected agate
Planes
Synthetic sapphire
(corundum)
Pans
Chromium plated brass

Pan supports
Fully adjustable
Arrestment
Drop action with front
release agate cones and
vees throughout
Case
Light grey enamel finish

Scale piece
40 divisions
Knife edges
Selected agate
Planes
Synthetic sapphire
(corundum)
Pans
Chromium plated brass

Pan supports
Fully adjustable
Arrestment
Drop action with front
release agate cones and
vees throughout
Case
Light grey enamel finish
with double front
counterpoised slides